

FULLY ADJUSTABLE INERT NEBULISER

FOR AA500 ATOMIC ABSORPTION SYSTEM

The NEW nebuliser system for the AA500 is available for retrofit to all existing instruments and will be supplied as standard on ALL new systems commencing October 2013.


This high precision nebuliser maximises stability and sensitivity even in the most complex matrices.

It is fully adjustable so that a wide variety of sample matrices – aqueous or organics (oils etc.), acid or alkali dilute or concentrated solutions can be analysed under optimum conditions.

- Fully inert nebuliser
- Adjustable from 2-6ml/min
- Excellent RSD's typically <1.0%
- Enhanced sensitivity
- Can be used on AA500 and AA990 systems
- Acid resistant version available with Pt/Ir capillary

We are offering a special introductory price for all retrofit nebuliser kits for the AA500 and AA990 until 31st December 2013.

Our retrofit update kits include nebuliser, modified spray chamber, fittings and instruction booklet.



PG Instruments Ltd, Alma Park, Wibtoft, Leicestershire, LE17 5BH, England

Tel: 0044 (0) 1455 220131

Fax: 0044 (0) 1455 220025

www.pg instruments.com

Email: sales@pginstruments.com

V1-10/13